# Quick Guide

# Spectroquant® Prove plus

Spectroquant® Prove

Spectrophotometer 100 plus • 300 plus • 600 plus

# Contents

1

<b>1</b>	<b>Safety .....</b>	<b>3</b>
	1.1 Intended Use .....	3
<b>2</b>	<b>Overview .....</b>	<b>4</b>
	2.1 Scope of Delivery .....	4
	2.2 Overview of the Instrument.....	4
	2.3 Display and User Interface .....	5
<b>3</b>	<b>Getting started .....</b>	<b>11</b>
	3.1 General Notes on Handling .....	11
	3.2 Initial Setup .....	11
	3.2.1 Connecting the Power Supply .....	11
	3.2.2 First Power-on .....	12
	3.2.3 Language Setup .....	13
	3.2.4 Date, Time and Country-specific Settings.....	13
	3.2.5 Self-test.....	14
<b>4</b>	<b>Operation.....</b>	<b>15</b>
	4.1 Switching the Spectrophotometer on or off.....	15
	4.2 System Setup .....	17
	4.3 Zero Adjustment for Preprogrammed Methods .....	18
	4.3.1 Notes on Zero Adjustment .....	19
	4.3.2 When to repeat the Zero Adjustment? .....	19
	4.4 Measurements.....	20
	4.4.1 Performing a Measurement .....	21
<b>5</b>	<b>Appendix .....</b>	<b>23</b>
	5.1 Certificates .....	23
	5.2 Warranty .....	24
	5.3 Accessories .....	25
	5.4 Table of Hazardous Substances' Name and Concentration ( <b>only for China</b> ).....	25

2

3

4

For more information about the Spectroquant® Prove plus  
please visit:  
[www.sigmaaldrich.com/spectroquant](http://www.sigmaaldrich.com/spectroquant)

5

For our Operating Manual and more information about  
the technical videos please visit:  
[www.sigmaaldrich.com/spectroquant](http://www.sigmaaldrich.com/spectroquant)

# 1 Safety

This Quick Guide contains basic instructions that you must follow during the commissioning, operation and maintenance of the spectrophotometer. Consequently, all responsible personnel must read this Quick Guide carefully before working with the meter. Keep this Quick Guide in the vicinity of the meter.

This is a class A device. This equipment may cause interference in a residential installation. In this case the user is encouraged to perform appropriate measures to correct the interference.

## Safety instructions

Safety instructions in this operating manual are indicated by the warning symbol (triangle) in the left margin. The signal word (such as "CAUTION") indicates the danger level.

The following warning symbols are used:

Symbols	Description
	<b>WARNING</b> Hazardous area (general). The xenon lamp (UV/VIS) emits radiation in the ultraviolet region, which may cause damage to the eyes. Never look directly in the radiation of this light source without wearing proper eye protection. Protect your skin from the direct exposure to UV light.
	<b>WARNING</b> Dangerous electrical voltage.
 <b>WARNING</b>	<b>WARNING</b> Signifies instructions that must be followed precisely in order to prevent serious dangers to personnel.
<b>CAUTION</b>	<b>CAUTION</b> Signifies instructions that must be followed precisely in order to avoid minor injuries to personnel or damage to the instrument or the environment.
 <b>CAUTION</b>	<b>CAUTION</b> This is a cautionary notice with a warning symbol drawing your attention to the risk of (limited) harm to personnel.
<b>NOTE</b>	<b>NOTE</b> Signifies a notice drawing your attention to special characteristics.
	<b>REFERENCE</b> Used to indicate references to other documents.

Please pay attention to the separate safety instructions leaflet (part of delivery scope) and read it carefully.

## 1.1 Intended Use

The intended use of the spectrophotometer consists exclusively of the carrying out of photometric measurements according to this operating manual. Observe the technical specifications of the cells in the operating manual.

Any other use is considered to be unauthorized. The spectrophotometer was developed for performing water analyses in the laboratory.

## 2 Overview

1

### 2.1 Scope of Delivery

- Spectrophotometer
- Power adapter
- Power connectors (3 pieces)
- Dust cover
- Zero cell
- Quick Guide (A4 format)
- Safety instructions
- Certificate of final inspection

2


### 2.2 Overview of the Instrument Packaging

The spectrophotometer is shipped in protective transport packaging.

#### CAUTION

Retain the original packaging including the inner packaging to protect the instrument against hard knocks if it has to be transported. Please note that damage caused by improper transport voids all warranty claims.

### Front of the instrument


3

4

5

## 2.3 Display and User Interface

### Ports at the rear of the instrument

- 5 Socket for plug-in power supply unit
- 6 LAN port
- 7 USB Mini B port
- 8 USB-A ports


### NOTE

All connections comply with SELV.

### NOTE

The entire display is touch-sensitive. Make selections using a fingertip or special touch pen. Do not touch the display with sharp objects (e.g. the tip of a ballpoint pen).

- Do not place objects on the display, as doing so may scratch it
- Touch buttons, words or symbols to select them
- Scrollbars are provided to assist quick movement through long lists
- Touch the arrow in the scrollbar to scroll upwards or downwards through the list
- Following selection, the item is activated immediately
- Touching a main button outlines it in blue
- Selecting an item inverts it (with dark text being shown on a light background)
- Selecting a text inverts it (with dark text being shown on a light background), e.g. method-specific settings for concentration mode "Show Absorbance"
- "0" is OFF, "I" is ON – the active selection is displayed light grey with dark figure, in this case the Show Absorbance is ON



1

**Main menu navigation**

The main menu is always visible on the left: It consist of two pages with four smart icons each. To switch between the two pages push  at the bottom on the left.

"Methods" and "Results" are the most often used modes and they are at the top of the main menu navigation.


2

3

**NOTE**

The menu selected is always outlined in blue.

4

**NOTE**

Action buttons like "Start", "Save", "Print" give the following touch feedback:

-  Normal  
Remains static

Active fields are always shown in bright color.

Pressed fields invert the color as long as the chosen action is performed.


-  Disabled  
Draws 30 % of the normal state

Inactive, disabled fields show faint color.

5

**NOTE**

The main menus "Settings (Method Settings)", "Ad hoc", "AQA", "System (Instrument Settings)", "Login/Logout", "Timer" open up a submenu. Example "Settings":


To leave these, the submenu has to be closed by touching the main menu button again, in this case:

The main menu "Methods" comprises two main overview panels arranged as shown below: the Concentration Measurement Overview and the Method List Overview.

**Screen layout concentration measurement overview**

1

Screen layout method list overview

The screenshot displays a software interface titled "Methods" with a dark theme. The interface includes a left sidebar with menu items: "Methods" (selected), "Settings", "Ad hoc", and "AQA". The main area shows a table of methods with columns for ID, Name, Range, Unit, and Application. Two dropdown menus at the top are both set to "All". A toolbar on the right contains search, refresh, and add buttons. A scrollbar is visible on the right side of the table.

Callouts and their corresponding UI elements:

- 1: Menu title (Methods)
- 2: Dropdown box (closed) (Left dropdown menu)
- 3: Dropdown box (closed) (Right dropdown menu)
- 4: Selection buttons (Search, Refresh, Add icons)
- 5: Scrollbar (Right side of the table)
- 6: Name of method (Aluminium in the table row)
- 7: Focuses selected main menu button (Methods icon in sidebar)

ID	Name	Range	Unit	Application
208	Acid Capacity	0.40-8.00 mg/l	OH	101758
2517	ADMI10	10-500	ADMI	App
2518	ADMI50	2.0-100.0	ADMI	App
196	Aluminium	0.02-0.50 mg/l	Al	100594
43	Aluminium	0.020-1.20 mg/l	Al	114825
2520	Ammonia, free	0.000-3.65 mg/l	NH <sub>3</sub>	App

2

3

4

5

## Overview of main buttons

Buttons	Description
	<b>Method list</b> List of all methods, irrespective of mode
	<b>Settings</b> This button is used to activate method-specific settings (e.g. sample dilution, turbidity correction, zero adjustment, sample blank, reagent blank)
	<b>Ad hoc</b> For performing measurements (absorbance/transmission, spectrum, kinetics) Allows measurements to be performed without the need to create methods
	<b>Absorbance/Transmission Mode</b> Ad hoc submenu: perform absorbance or transmission measurements
	<b>Spectrum Mode</b> Ad hoc submenu: record spectrum Method list: create methods -> Spectrum Mode
	<b>Kinetic Mode</b> Ad hoc submenu: perform kinetic measurement Method list: create methods -> Kinetic Mode
	<b>AQA</b> Overview and list of all Analytical Quality Assurance (AQA) modes
	<b>AQA Status 1&amp;2</b> AQA submenus: Status display of the period of validity and the outcome (passed/failed)
	<b>AQA1</b> AQA submenu: List of AQA1 methods
	<b>AQA2</b> AQA submenu: List of AQA2 methods
	<b>Pipette check</b> AQA submenu: List of pipette-checking methods
	<b>Result list</b> List of all stored results
	<b>System setup</b> This button is for optional instrument settings (e.g. date, time, updates etc.)
	<b>Login/logout</b> Log users in and out
	<b>Timer list</b> List of stopwatch functions

1

Overview of action & selection buttons

Action & Selection Buttons	Description
	<b>Start button</b> Start an action (e.g. measurement)
	<b>Start zero</b> Start zero adjustment for a method
	<b>Apply</b>
	<b>Save</b>
	<b>Stop</b>
	<b>Close</b>
	<b>Logout</b> User logout
	<b>Search method</b>
	<b>Search/results list</b> Search function, search criterion: method name, method number or item number
	<b>Filter cancellation button</b> Cancel all set filter options
	<b>Edit</b> For editing parameters
	<b>Create method</b>
	<b>Print</b> Print to pdf (USB device) or printer
	<b>Export button</b> All selected results are exported to an external memory device as .csv file
	<b>Import button</b> Updates/Methods are imported from an external memory device into the instrument
	<b>Delete</b> The selected items are deleted

2

3

4

5

# 3 Getting Started

## 3.1 General Notes on Handling

The Spectroquant® Prove plus spectrophotometer is an optical precision instrument. It should always be handled with care, especially when in mobile use. Always protect the instrument from conditions that could damage the mechanical, optical and electrical components. Please note the following in particular:

- The temperature and humidity during operation and storage must be within the limits specified in the "Technical Data" section (see the Operating Manual)

**The instrument must never be exposed to the following:**

- Extreme dust, humidity and moisture
- Intense light and heat
- Fumes that are corrosive or contain high concentrations of solvents

**In addition take care of the following:**

- For measuring, the instrument must be placed on a flat surface
- Spilled liquid or other material should be removed immediately (see the Operating Manual)
- If a cell has broken in the cell holder, the cell holder should be cleaned immediately (see the Operating Manual)
- The cover should always be closed when the spectrophotometer is not in use
- When the spectrophotometer is being transported, the cell compartment must be empty

## 3.2 Initial Setup

**Proceed as follows:**

- Connect the power adapter (see section 3.2.1)
- Switch on the spectrophotometer (see section 3.2.2)
- Set the language (see section 3.2.3)
- Set the date and time (see section 3.2.4)
- Run the self-test (see section 3.2.5)

### NOTE

For our Operating Manual and more information about the technical videos please visit: [www.sigmaaldrich.com/spectroquant](http://www.sigmaaldrich.com/spectroquant)

### 3.2.1 Connecting the Power Supply

Power is supplied through the power adapter provided. The power adapter supplies the spectrophotometer with the required voltage and type of current (24 V DC).

### CAUTION

The line voltage at the user location must fulfil the specifications stated on the power adapter (the specifications are also indicated in the Operating Manual). Only ever use the 24 V power adapter provided. Please note that damage caused by using a different power adapter than the one supplied voids all warranty claims.

1


2

3

4

5


1


2

### 3.2.2 First Power-on

After switching on the spectrophotometer for the first time you are automatically guided through the language, date and time setup procedures.

3

1. Press the ON/OFF button **1**. The spectrophotometer gives an audible signal (beep) and starts booting for approximately 2 minutes. You will see the following display:

#### Connecting the power adapter:

1. Connect the miniplug **1** of the power adapter to the socket **2** of the spectrophotometer.
2. Connect the power adapter to a wall socket.


2. The display switches to language setup (see section 3.2.3).

4

5

### 3.2.3 Language Setup

The software supports several languages. When you switch on the spectrophotometer for the first time, a list of language options is automatically displayed after the boot procedure.


1. Select the desired language (1).
2. Tap on the Save button (2) to confirm.

#### NOTE


The saving process of changing the language requires some seconds.

### 3.2.4 Date, Time and Country-specific Settings

During initial setup, having set the language option you are automatically guided through the date and time setup procedure.


1. Tap on the Date format button (1).
2. The calendar view pops up (2). You can now enter the date.

3. Tap on OK (3) to confirm.
4. Tap on the Arrow button (4) to choose the country-specific basic date setting. The date format can be set and displayed for EU and US.

1

2

5. Tap on the Time format button **5**. The numeric key panel **6** pops up. Now you can enter the time.
6. Tap on OK **7** to confirm.
7. Tap on the Arrow button **8** to choose the country-specific basic time setting. The time format can be set and displayed for EU and US.
8. Tap on the Arrow button **9** to choose the decimal separator "."/"/," used in your country.
9. Tap on the Save button **10** to confirm.


3

4

5

### 3.2.5 Self-test

Following language, date and time setup the spectrophotometer performs a self-test.


1. Remove all cells and close the cell compartment cover.
2. Start the self-test with the Start button **1**.
3. The spectrophotometer performs the self-test.

#### Self-test

The self-test covers:


- Checks on memory, processor, internal interfaces, filter and lamp
- A calibration of the wavelength

When the self-test has ended, the display shows the main menu.


# 4 Operation

## 4.1 Switching the Spectrophotometer On or Off

### Switching on


1. Press the ON/OFF button **1**. The spectrophotometer gives an audible signal (beep) and starts booting for approximately 2 minutes. You will see the following display:

2. After the booting process the screen shows the self-test dialog.

### Starting the self-test

3. Remove all cells and close the cell compartment cover.

4. Start the self-test with the Start button **2**.
5. The spectrophotometer performs the self-test.

### Self-test

The self-test covers:


- Checks on memory, processor, internal interfaces, filter and lamp
- A calibration of the wavelength

When the self-test has ended, the display shows the main menu.

1


### Energy-saving mode – display

2

### Switching off

Press the ON/OFF button 4 to switch the spectrophotometer off.

3

The spectrophotometer automatically switches off the backlight of the display 3 when no button has been tapped within a period of 10 minutes. The backlight is switched on again with the next tap. The button functions are activated only following a further tap.

### NOTE

The instrument has an Auto-Power-Off function, which switches it automatically off after a user-defined time. This function is not active out of the box, but you can turn it on in "System (Instrument settings)".

### NOTE

You can set a user-defined time for this function (please see the Operating Manual for detailed instructions).

4


5

## 4.2 System Setup

General instrument setup is carried out in the "System" menu.

Buttons	Description
	<p><b>Information</b></p> <p>This submenu displays the following information about the device: Software/method versions, device class, lamp counter and serial number</p>
	<p><b>Interface</b></p> <p>This submenu displays the following settings options and standard settings: Audible signals – ON, Backlight – 100 %, Print to pdf – ON</p>
	<p><b>Region</b></p> <p>This submenu displays the following settings options and standard settings: Language, date, time and country zone EU/US, decimal separator – "." / "," (dot or comma)</p>
	<p><b>Quality</b></p> <p>This submenu displays the following settings options and standard settings: Quick zero – OFF, AQA1 and AQA2 lock – OFF, Zero Adjustment expiry – ON (interval: 7 days), Use expired reagents – OFF, Service reminder – ON</p>
	<p><b>Automation (Setup 4)</b></p> <p>This submenu displays the following settings options and standard settings: Energy saving mode – ON (10 minutes), Auto Power off – OFF, Auto log off – OFF, Auto store – ON, Auto print – OFF, Sample ID popup – OFF</p>
	<p><b>User management</b></p> <p>This submenu displays the following settings options and standard settings: Activation of user management and administrator settings, User login required – OFF</p>
	<p><b>Service</b></p> <p>This submenu displays the following settings options: Various service functions such as backup, restore, export of log or system data and import of methods</p>
	<p><b>Update</b></p> <p>This submenu displays the option for performing software and method updates</p>

1

Buttons	Description
	<b>Network</b> This submenu displays the setting options for connecting the Prove plus device with a network
	<b>Prove Connect</b> This submenu displays the settings options for connecting the Prove plus device with the Prove Connect software (the Prove Connect software is optionally available, order No. Prove Connect to LIMS Y110860001)

2

Please see the Operating Manual for detailed information on the buttons and their functions.


3

### 4.3 Zero Adjustment for Preprogrammed Methods


A zero adjustment must be performed for each cell type. The zero adjustment for concentration methods is stored within the spectrophotometer separately for each cell type. The period of validity of the zero adjustment for concentration methods can be edited in the System settings. When a zero adjustment has already been performed for the inserted cell type and the selected method, the date of the most recent zero adjustment is displayed in the info line.

When a zero adjustment is required, following screen will pop up:

4

5

1. Insert the zero cell according to cell type. Zero adjustment starts automatically and, if the zero adjustment is passed, a tick **1** appears in the status display field for the Zero adjustment.  
In the case of a method which only measures the sample at a single wavelength the absorbance of the Zero value **2** is also displayed.
2. With a cell inserted, the zero adjustment can be repeated manually by tapping on the Start zero button **3**.
3. Tapping on the OK button accepts the zero adjustment value for the method.
4. The screen changes to show the concentration measurement screen (see page 7).
5. The instrument is ready to start measuring the sample.

**NOTE**

The instrument provides a Quick Zero function. With this the zero adjustment is performed for all wavelengths which are used with all Spectroquant® test kits. This function can be activated via "System (Instrument settings)" – "Quality".

**4.3.1 Notes on Zero Adjustment****Zero adjustment with round cells**

- Only use clean, scratch-free round cells and distilled water.  
The minimum filling level is 20 mm. A ready zero cell is contained within the scope of delivery of the spectrophotometer
- A ready zero cell can, in principle, be used for an indefinite period of time. We recommend, however, that you regularly check the zero cell for visible contamination and scratches and refill or exchange it if necessary (at least every 24 months)
- Insert the round cell until it touches the bottom of the round cell compartment

**Zero adjustment with rectangular cells**

- With rectangular cells, zero adjustment must be carried out using the same cell type (manufacturer and cell material [e.g. optical glass, quartz glass, plastic]) as the one that will be used for measurement. This is important because cells of different manufacturers have different absorption characteristics. When changing the cell type, repeat the zero adjustment with the new type
- Prior to zero adjustment, clean the rectangular cell and fill it with distilled water. The minimum filling level is 20 mm
- Rectangular cells always have to be inserted in the cell compartment with the same orientation for measurement and zero adjustment (e.g. cell inscription always on the left side)
- Insert the rectangular cell until it touches the bottom and left edge of the holder. The opaque sides of the rectangular cell must point to the front and rear

**4.3.2 When to repeat the Zero Adjustment?****We recommend that you repeat the zero adjustment in the following cases:**

- If the spectrophotometer was subject to mechanical stress such as strong shock or transport
- If the ambient temperature has changed by more than 5 °C since the last zero adjustment
- At least once a week. The interval to repeat a zero adjustment is set in the instrument to 7 days. You can change this under "System (Instrument Settings)"
- If a new cell type (different manufacturer, different glass type) is used
- Basically, each time you want to measure with the highest possible accuracy

**NOTE**

If an interval to repeat a zero adjustment is set you will be prompted to repeat it after the interval has passed. You can also repeat a zero adjustment by selecting a method, then touching the "Settings" icon. Choose "Zero adjustment" and insert a zero cell to start the measurement.

1

2

3

4

5

1

## 4.4 Measurements

The spectrophotometer can be used to perform the measurements listed below.

2

Type of Measurement	Description
Concentration	<ul style="list-style-type: none"> <li>• Preprogrammed methods that can be executed using Spectroquant® test kits or self-prepared reagents</li> <li>• User-programmed methods</li> </ul>
Absorption / transmission	<ul style="list-style-type: none"> <li>• Single-wavelength measurements for establishing the absorbance or transmission of solutions</li> <li>• Multiple-wavelength measurements for establishing the absorbance or transmission of solutions</li> </ul>
Spectrum	<ul style="list-style-type: none"> <li>• Programmed methods for establishing the absorbance or transmission of solutions over a defined wavelength range</li> </ul>
Kinetics	<ul style="list-style-type: none"> <li>• Programmed methods for establishing the absorbance or transmission of solutions over a defined period</li> </ul>
Quality checks	Instrument supported analytical quality assurance: <ul style="list-style-type: none"> <li>• Instrument check (AQA 1)</li> <li>• Method-specific system check – preprogrammed for all Spectroquant® standards (AQA 2)</li> <li>• Pipette volume control (PipeCheck)</li> <li>• Check interferences from foreign substances (MatrixCheck)</li> </ul>

3


4

5

### 4.4.1 Performing a Measurement


Measurements can be performed using rectangular cells of various path lengths (10, 20, 50 mm/100 mm Prove 600 plus) and Spectroquant® round cells. Insert cells as follows to start the measurement:

#### Measuring with a round cell with closed lid

- Insert the barcoded Spectroquant® round cell through the opening **1**, ensuring that the white position mark **2** on the cell is aligned with the positioning mark on the spectrophotometer **3**
- Measurement starts automatically, and the measurement result is displayed in the concentration measurement overview (see page 7)

#### Measuring with rectangular cells with open lid: insert AutoSelector


- Open the flip-up cover **4** by pushing it back with your fingers
- Insert the AutoSelector vertically into the cell compartment **5**, ensuring that the white position mark **6** on the AutoSelector is aligned with the positioning mark on the spectrophotometer **7**
- The photometer is ready to measure

#### NOTE

If the barcode cannot be read, please see the [Operating Manual](#).

1


**Measuring with rectangular cells with open lid:  
Insert rectangular cells (10, 20, 50 mm)**

- Insert the rectangular cell vertically into the cell compartment **8**, ensuring that the cell is flush against the left side of the cell holder **9** at all times
- Measurement starts automatically, and the measurement result is displayed in the concentration measurement overview (see page 7)

2

**Measuring with rectangular cells with open lid:  
Insert 100 mm rectangular cells (Prove 600 plus)**

- Remove the top of the round cell compartment including the AutoSelector **10**
- Insert the 100 mm rectangular cell vertically into the cell holder **11**. Make sure that you hold it with both hands on the small edges while inserting it carefully
- Measurement starts automatically, and the measurement result is displayed in the concentration measurement overview (see page 7)

**NOTE**

Please see Analytical Procedures and Appendices for detailed measurement procedures.

3

4

5

# 5 Appendix

## 5.1 Certificates

		Ref. Certif. No. <b>DE 3 - 32051</b>
IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE) CB SCHEME		
<b>CB TEST CERTIFICATE</b>		
Product	Laboratory equipment	
Name and address of the applicant	<b>Merck KGaA</b> Frankfurter Straße 250 64293 Darmstadt GERMANY	
Name and address of the manufacturer	<b>Merck KGaA</b> Frankfurter Straße 250, 64293 Darmstadt, GERMANY	
Name and address of the factory	<b>LRE Medical GmbH</b> Hofer Straße 5, 86720 Nördlingen, GERMANY	
Ratings and principal characteristics	Rated voltage: 24 VDC Rated frequency: 2.5 A	
Trademark / Brand (if any)	Merck	
Model/type Ref.	Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600	
Additional information (if necessary)	Certificate DE 3 - 3922 issued 2014-11-13 is replaced by this version due to technical changes.	
A sample of the product was tested and found to be in conformity with	IEC 61010-1:2010 IEC 61010-1:2010/AMD1:2016	
as shown in the Test Report Ref. No. which forms part of this certificate	713298142	
Page 1 of 1 This CB Test Certificate is issued by the National Certification Body		
CB 034009 0010 Rev. 01 Date: 2023-05-05		
		
TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany		

Rev. 00		Test Report issued under the responsibility of: <b>NCB TÜV SÜD Product Service GmbH</b> Ridlerstraße 65, D-80339 München Germany	
<b>TEST REPORT</b> <b>IEC 61010-1</b> <b>Safety requirements for electrical equipment for measurement, control, and laboratory use</b> <b>Part 1: General requirements</b>			
Report Number.	713298142		
Date of issue	2023-04-28		
Total number of pages	79		
Name of Testing Laboratory preparing the Report	TÜV SÜD Product Service GmbH		
Applicant's name	Merck KGaA		
Address	Frankfurter Str. 250, DE-64293 Darmstadt, Germany		
Test specification:	Standard ..... IEC 61010-1:2010, IEC 61010-1:2010/AMD1:2016 Test procedure..... CB Scheme Non-standard test method ..... N/A		
TRF template used	IECEE OD-2020-F1:2020, Ed.1.3		
Test Report Form No.	IEC61010_1P		
Test Report Form(s) Originator	VDE Prüf- und Zertifizierungsinstitut GmbH		
Master TRF	2021-04-12		
Copyright © 2021 IEC System of Conformity Assessment Schemes for Electrotechnical Equipment and Components (IECEE System). All rights reserved. This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context. If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo and the reference to the CB Scheme procedure shall be removed. <b>This report is not valid as a CB Test Report unless signed by an approved IECEE Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.</b>			
<b>General disclaimer:</b> The test results presented in this report relate only to the object tested. This report shall not be reproduced, except in full, without the written approval of the Issuing NCB. The authenticity of this Test Report and its contents can be verified by contacting the NCB, responsible for this Test Report.			

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證證書 ◆ СЕРТИФИКАТ ◆ CERTIFICADO ◆ CERTIFICAT			
	<b>CERTIFICATE</b> No. U8 034909 0011 Rev. 01		
	Holder of Certificate:	<b>Merck KGaA</b> Frankfurter Straße 250 64293 Darmstadt GERMANY	
	Certification Mark:		
	Product:	Laboratory Equipment	
	Brand Name:	Merck	
	Model(s):	Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600	
	Parameters:	Rated voltage: 24 VDC Rated frequency: 2.5 A	
	Tested according to:	UL 61010-1:2012/R:2019-07 CSA C22.2 No. 61010-1:2012/A1:2018-11	
	This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL for USA and a Standards Council of Canada ISO/IEC 17065 accredited Certification body for Canada. Test report no.: 713298142		
Date:	2023-05-08		
			
TÜV SÜD America, Inc. • 401 Edgewater Place Suite #600 • Wakefield • MA 01880 • USA			

<b>Technical Report</b>		
Technical Report No.: 713298142 Date: 2023-04-28		Product Service
Client:	LRE Medical GmbH  Georg-Brauchle-Ring 89, DE-80992 München Germany	
Licence holder:	Merck KGaA  Frankfurter Str. 250, DE-64293 Darmstadt Germany	
Factory:	LRE Medical GmbH  Hofer Straße 5, DE-86720 Nördlingen Germany	
Test object:	Product: Laboratory equipment Model: Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600	
Test specification:	IEC 61010-1:2010 IEC 61010-1:2010/AMD1:2016  EN 61010-1:2010/A1:2019  UL 61010-1:2012/R:2019-07 CSA C22.2 No. 61010-1:2012/A1:2018-11	
Purpose of examination:	Testing and evaluation according to the test specification	
Test result:	The test results show that the presented product is in compliance with the above listed test specifications.	
Any use for advertising purposes must be granted in writing. This technical report may only be quoted in full. This report is the result of a single examination of the object in question. It does not imply a general statement regarding the quality of products from regular production. For further details please see testing and certification regulation, chapter A-3.4.		
Report No.: 713298142	Project Manager: Ralph Fischer	Phone: +49 89 5098-4322
Rev. 00	Date: 2023-04-28	Fax: +49 89 5098-4130
Page 1 of 4		E-Mail: ralph.fischer@tvsud.com
		30359 München
		Deutschland

1

## 5.2 Warranty

In addition to any other rights provided by local law which shall not be limited herewith, the manufacturer provides an additional limited worldwide warranty for a period of 12 months from the date of shipment, that the spectrophotometer is free from any defects that are due to faulty material or workmanship unless otherwise noted in the product manual. This warranty does not apply to consumable products such as the halogen lamp.

2

The spectrophotometer must only be opened, adjusted or repaired by specialist personnel authorized by the manufacturer. Noncompliance invalidates any warranty claims.

Keep the original packing including the inner packing to protect the instrument against hard shocks if it has to be shipped. Note that damage caused by improper transport voids all warranty claims.

3

In case of a warranty claim please contact the technical customer service of your local supplier to get further information about the required documents and possibilities of warranty service (e.g. replacement, repair) if a warranty claim is approved.

4

5

### 5.3 Accessories

Description	Order No.
Halogen lamp module for Spectroquant® Prove 100	1.74010.0001
Case for Spectroquant® spectrophotometer Prove 100   300 and 600	1.73020.0001
Rectangular cells 10 mm (1 pack = 2 pcs)	1.14946.0001
Rectangular cells 20 mm (1 pack = 2 pcs)	1.14947.0001
Rectangular cells 50 mm (1 pack = 2 pcs)	1.14944.0001
Semi-microcells 50 mm (1 pack = 2 pcs)	1.73502.0001
Rectangular cells quartz 10 mm (1 pack = 2 pcs)	1.00784.0001
Empty cells 16 mm Ø (1 pack = 25 pcs) with screw cap	1.14724.0001
Zero Cell (1 pack = 1 pc)	1.73503.0001
Rectangular cell 100 mm	1.74011.0001
Prove Connect to LIMS Unlimited License	Y110860001

### 5.4 Table of Hazardous Substances' Name and Concentration (only for China)

Spectroquant® Prove 100 plus | 300 plus | 600 plus systems

Component Name	Hazardous substances' name					
	(Pb)	(Hg)	(Cd)	(Cr <sup>6+</sup> )	(PBB)	(PBDE)
Housing	O	O	O	O	O	O
Optics	O	O	O	O	O	O
Electronics	O	O	O	O	O	O

O: Indicates that this toxic or hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in SJ / T11363-2006.

X: Indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in SJ / T11363-2006.

- Data listed in the table represents best information available at the time of publication
- Product to be operated under Environmental Conditions as defined in the Product Manual to maintain the declared EFUP
- Consumables or parts may have their own label with an EFUP value less than the system. Periodic replacement of those consumables or parts to maintain the declared EFUP is done in accordance with the declared EFUP is done in accordance with the Product Maintenance Procedures

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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# 5 Anhang

## 5.1 Zertifikate

		Ref. Certif. No. <b>DE 3 - 32051</b>
IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE) CB SCHEME		
<b>CB TEST CERTIFICATE</b>		
Product	Laboratory equipment	
Name and address of the applicant	Merck KGaA Frankfurter Straße 250 64293 Darmstadt GERMANY	
Name and address of the manufacturer	Merck KGaA Frankfurter Straße 250, 64293 Darmstadt, GERMANY	
Name and address of the factory	LRE Medical GmbH Hofer Straße 5, 86720 Nördlingen, GERMANY	
Ratings and principal characteristics	Rated voltage: 24 VDC Rated frequency: 2.5 A	
Trademark / Brand (if any)	Merck	
Model/type Ref.	Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600	
Additional information (if necessary)	Certificate DE 3 - 3922 issued 2014-11-13 is replaced by this version due to technical changes.	
A sample of the product was tested and found to be in conformity with	IEC 61010-1:2010 IEC 61010-1:2010/AMD1:2016	
as shown in the Test Report Ref. No. which forms part of this certificate	713298142	
Page 1 of 1 This CB Test Certificate is issued by the National Certification Body		
CB 034009 0010 Rev. 01 Date: 2023-05-05		
		
TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany		

Rev. 00		Test Report issued under the responsibility of: NCB TÜV SÜD Product Service GmbH Ridlerstraße 65, D-80339 München Germany	
<b>TEST REPORT</b> <b>IEC 61010-1</b> <b>Safety requirements for electrical equipment for measurement, control, and laboratory use</b> <b>Part 1: General requirements</b>			
Report Number.	713298142		
Date of issue	2023-04-28		
Total number of pages	79		
Name of Testing Laboratory preparing the Report	TÜV SÜD Product Service GmbH		
Applicant's name	Merck KGaA		
Address	Frankfurter Str. 250, DE-64293 Darmstadt, Germany		
Test specification:	Standard ..... IEC 61010-1:2010, IEC 61010-1:2010/AMD1:2016 Test procedure..... CB Scheme Non-standard test method ..... N/A		
TRF template used	IECEE OD-2020-F1:2020, Ed.1.3		
Test Report Form No.	IEC61010_1P		
Test Report Form(s) Originator	VDE Prüf- und Zertifizierungsinstitut GmbH		
Master TRF	2021-04-12		
Copyright © 2021 IEC System of Conformity Assessment Schemes for Electrotechnical Equipment and Components (IECEE System). All rights reserved. This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context. If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo and the reference to the CB Scheme procedure shall be removed. <b>This report is not valid as a CB Test Report unless signed by an approved IECEE Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.</b>			
<b>General disclaimer:</b> The test results presented in this report relate only to the object tested. This report shall not be reproduced, except in full, without the written approval of the Issuing NCB. The authenticity of this Test Report and its contents can be verified by contacting the NCB, responsible for this Test Report.			

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證證書 ◆ СЕРТИФИКАТ ◆ CERTIFICADO ◆ CERTIFICAT			
	<b>CERTIFICATE</b> No. U8 034909 0011 Rev. 01		
	Holder of Certificate:	Merck KGaA Frankfurter Straße 250 64293 Darmstadt GERMANY	
	Certification Mark:		
	Product:	Laboratory Equipment	
	Brand Name:	Merck	
	Model(s):	Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600	
	Parameters:	Rated voltage: 24 VDC Rated frequency: 2.5 A	
	Tested according to:	UL 61010-1:2012/R:2019-07 CSA C22.2 No. 61010-1:2012/A1:2018-11	
	This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL for USA and a Standards Council of Canada ISO/IEC 17065 accredited Certification body for Canada. Test report no.: 713298142		
Date:	2023-05-08		
			
TÜV SÜD America, Inc. • 401 Edgewater Place Suite #600 • Wakefield • MA 01880 • USA			

<b>Technical Report</b>		
Technical Report No.: 713298142 Date: 2023-04-28		Product Service
Client:	LRE Medical GmbH  Georg-Brauchle-Ring 89, DE-80992 München Germany	
Licence holder:	Merck KGaA  Frankfurter Str. 250, DE-64293 Darmstadt Germany	
Factory:	LRE Medical GmbH  Hofer Straße 5, DE-86720 Nördlingen Germany	
Test object:	Product: Laboratory equipment Model: Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600	
Test specification:	IEC 61010-1:2010 IEC 61010-1:2010/AMD1:2016  EN 61010-1:2010/A1:2019  UL 61010-1:2012/R:2019-07 CSA C22.2 No. 61010-1:2012/A1:2018-11	
Purpose of examination:	Testing and evaluation according to the test specification	
Test result:	The test results show that the presented product is in compliance with the above listed test specifications.	
Any use for advertising purposes must be granted in writing. This technical report may only be quoted in full. This report is the result of a single examination of the object in question. It does not imply a general statement regarding the quality of products from regular production. For further details please see testing and certification regulation, chapter A-3.4.		
Report No.: 713298142	Project Manager: Ralph Fischer	Phone: +49 89 5098-4322
Rev: 00	Date: 2023-04-28	Fax: +49 89 5098-4130
Page 1 of 4		E-Mail: ralph.fischer@tvsud.com
		80339 München
		Deutschland

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# 5 Annexes

## 5.1 Certificats

		Ref. Certif. No. <b>DE 3 - 32051</b>
IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE) CB SCHEME		
<b>CB TEST CERTIFICATE</b>		
Product	Laboratory equipment	
Name and address of the applicant	<b>Merck KGaA</b> Frankfurter Straße 250 64293 Darmstadt GERMANY	
Name and address of the manufacturer	Merck KGaA Frankfurter Straße 250, 64293 Darmstadt, GERMANY	
Name and address of the factory	LRE Medical GmbH Hofer Straße 5, 86720 Nördlingen, GERMANY	
Ratings and principal characteristics	Rated voltage: 24 VDC Rated frequency: 2.5 A	
Trademark / Brand (if any)	Merck	
Model/type Ref.	Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600	
Additional information (if necessary)	Certificate DE 3 - 3922 issued 2014-11-13 is replaced by this version due to technical changes.	
A sample of the product was tested and found to be in conformity with	IEC 61010-1:2010 IEC 61010-1:2010/AMD1:2016	
as shown in the Test Report Ref. No. which forms part of this certificate	713298142	
Page 1 of 1 This CB Test Certificate is issued by the National Certification Body		
CB 034009 0010 Rev. 01 Date: 2023-05-05		
		
TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany		

Rev. 00 	Test Report issued under the responsibility of: <b>NCB TÜV SÜD Product Service GmbH</b> Ridlerstraße 65, D-80339 München Germany	
<b>TEST REPORT</b> <b>IEC 61010-1</b> <b>Safety requirements for electrical equipment for measurement, control, and laboratory use</b> <b>Part 1: General requirements</b>		
Report Number. ....	713298142	
Date of issue .....	2023-04-28	
Total number of pages.....	79	
Name of Testing Laboratory preparing the Report.....	TÜV SÜD Product Service GmbH	
Applicant's name.....	Merck KGaA	
Address.....	Frankfurter Str. 250, DE-64293 Darmstadt, Germany	
Test specification:		
Standard .....	IEC 61010-1:2010, IEC 61010-1:2010/AMD1:2016	
Test procedure.....	CB Scheme	
Non-standard test method .....	N/A	
TRF template used .....	IECEE OD-2020-F1:2020, Ed.1.3	
Test Report Form No. ....	IEC61010_1P	
Test Report Form(s) Originator .....	VDE Prüf- und Zertifizierungsinstitut GmbH	
Master TRF .....	2021-04-12	
Copyright © 2021 IEC System of Conformity Assessment Schemes for Electrotechnical Equipment and Components (IECEE System). All rights reserved. This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context. If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo and the reference to the CB Scheme procedure shall be removed. <b>This report is not valid as a CB Test Report unless signed by an approved IECEE Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.</b>		
<b>General disclaimer:</b> The test results presented in this report relate only to the object tested. This report shall not be reproduced, except in full, without the written approval of the Issuing NCB. The authenticity of this Test Report and its contents can be verified by contacting the NCB, responsible for this Test Report.		

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證證書 ◆ CERTIFICADO ◆ CERTIFICAT ◆ CERTIFIKAT ◆ CERTIFICATE				
	<b>CERTIFICATE</b> No. U8 034909 0011 Rev. 01			
	Holder of Certificate:	<b>Merck KGaA</b> Frankfurter Straße 250 64293 Darmstadt GERMANY		
	Certification Mark:			
	Product:	Laboratory Equipment		
	Brand Name:	Merck		
	Model(s):	Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600		
	Parameters:	Rated voltage: 24 VDC Rated frequency: 2.5 A		
	Tested according to:	UL 61010-1:2012/R:2019-07 CSA C22.2 No. 61010-1:2012/A1:2018-11		
	This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL for USA and a Standards Council of Canada ISO/IEC 17065 accredited Certification body for Canada. Test report no.: 713298142			
Date:	2023-05-08			
				
TÜV SÜD America, Inc. • 401 Edgewater Place Suite #600 • Wakefield • MA 01880 • USA				

<b>Technical Report</b>		
Technical Report No.: 713298142 Date: 2023-04-28		
Client:	LRE Medical GmbH  Georg-Brauchle-Ring 89, DE-80992 München Germany	
Licence holder:	Merck KGaA  Frankfurter Str. 250, DE-64293 Darmstadt Germany	
Factory:	LRE Medical GmbH  Hofer Straße 5, DE-86720 Nördlingen Germany	
Test object:	Product: Laboratory equipment Model: Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600	
Test specification:	IEC 61010-1:2010 IEC 61010-1:2010/AMD1:2016  EN 61010-1:2010/A1:2019  UL 61010-1:2012/R:2019-07 CSA C22.2 No. 61010-1:2012/A1:2018-11	
Purpose of examination:	Testing and evaluation according to the test specification	
Test result:	The test results show that the presented product is in compliance with the above listed test specifications.	
Any use for advertising purposes must be granted in writing. This technical report may only be quoted in full. This report is the result of a single examination of the object in question. It does not imply a general statement regarding the quality of products from regular production. For further details please see testing and certification regulation: chapter A-3.4.		
Report No.: 713298142	Project Manager: Ralph Fischer	Phone: +49 89 5098-4322
Rev: 00	Date: 2023-04-28	Fax: +49 89 5098-4130
Page 1 of 4		E-Mail: ralph.fischer@tvsud.com
		30359 München
		Deutschland

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Rev. 00 	Test Report issued under the responsibility of: NCB TÜV SÜD Product Service GmbH Ridlerstraße 65, D-80339 München Germany	
<b>TEST REPORT</b> <b>IEC 61010-1</b> <b>Safety requirements for electrical equipment for measurement, control, and laboratory use</b> <b>Part 1: General requirements</b>		
Report Number. ....	713298142	
Date of issue .....	2023-04-28	
Total number of pages.....	79	
Name of Testing Laboratory preparing the Report.....	TÜV SÜD Product Service GmbH	
Applicant's name.....	Merck KGaA	
Address.....	Frankfurter Str. 250, DE-64293 Darmstadt, Germany	
Test specification:		
Standard .....	IEC 61010-1:2010, IEC 61010-1:2010/AMD1:2016	
Test procedure.....	CB Scheme	
Non-standard test method .....	N/A	
TRF template used .....	IECEE OD-2020-F1:2020, Ed.1.3	
Test Report Form No. ....	IEC61010_1P	
Test Report Form(s) Originator .....	VDE Prüf- und Zertifizierungsinstitut GmbH	
Master TRF .....	2021-04-12	
Copyright © 2021 IEC System of Conformity Assessment Schemes for Electrotechnical Equipment and Components (IECEE System). All rights reserved. This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context. If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo and the reference to the CB Scheme procedure shall be removed. <b>This report is not valid as a CB Test Report unless signed by an approved IECEE Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.</b>		
<b>General disclaimer:</b> The test results presented in this report relate only to the object tested. This report shall not be reproduced, except in full, without the written approval of the Issuing NCB. The authenticity of this Test Report and its contents can be verified by contacting the NCB, responsible for this Test Report.		

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證證書 ◆ СЕРТИФИКАТ ◆ CERTIFICADO ◆ CERTIFICAT				
	<b>CERTIFICATE</b> No. U8 034909 0011 Rev. 01			
	Holder of Certificate:	Merck KGaA Frankfurter Straße 250 64293 Darmstadt GERMANY		
	Certification Mark:			
	Product:	Laboratory Equipment		
	Brand Name:	Merck		
	Model(s):	Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600		
	Parameters:	Rated voltage: 24 VDC Rated frequency: 2.5 A		
	Tested according to:	UL 61010-1:2012/R:2019-07 CSA C22.2 No. 61010-1:2012/A1:2018-11		
	This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL for USA and a Standards Council of Canada ISO/IEC 17065 accredited Certification body for Canada. Test report no.: 713298142			
Date:	2023-05-08			
				
TÜV SÜD America, Inc. • 401 Edgewater Place Suite #600 • Wakefield • MA 01880 • USA				

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Rev. 00 	Test Report issued under the responsibility of: NCB TÜV SÜD Product Service GmbH Ridlerstraße 65, D-80339 München Germany	
<b>TEST REPORT</b> <b>IEC 61010-1</b> <b>Safety requirements for electrical equipment for measurement, control, and laboratory use</b> <b>Part 1: General requirements</b>		
Report Number.	713298142	
Date of issue	2023-04-28	
Total number of pages	79	
Name of Testing Laboratory preparing the Report	TÜV SÜD Product Service GmbH	
Applicant's name	Merck KGaA	
Address	Frankfurter Str. 250, DE-64293 Darmstadt, Germany	
Test specification:		
Standard	IEC 61010-1:2010, IEC 61010-1:2010/AMD1:2016	
Test procedure	CB Scheme	
Non-standard test method	N/A	
TRF template used	IECEE OD-2020-F1:2020, Ed.1.3	
Test Report Form No.	IEC61010_1P	
Test Report Form(s) Originator	VDE Prüf- und Zertifizierungsinstitut GmbH	
Master TRF	2021-04-12	
Copyright © 2021 IEC System of Conformity Assessment Schemes for Electrotechnical Equipment and Components (IECEE System). All rights reserved. This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context. If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo and the reference to the CB Scheme procedure shall be removed. <b>This report is not valid as a CB Test Report unless signed by an approved IECEE Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.</b>		
<b>General disclaimer:</b> The test results presented in this report relate only to the object tested. This report shall not be reproduced, except in full, without the written approval of the Issuing NCB. The authenticity of this Test Report and its contents can be verified by contacting the NCB, responsible for this Test Report.		

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證證書 ◆ СЕРТИФИКАТ ◆ CERTIFICADO ◆ CERTIFICAT				
	<b>CERTIFICATE</b> No. U8 034909 0011 Rev. 01			
	Holder of Certificate:	Merck KGaA Frankfurter Straße 250 64293 Darmstadt GERMANY		
	Certification Mark:			
	Product:	Laboratory Equipment		
	Brand Name:	Merck		
	Model(s):	Spectroquant Prove 100, Spectroquant Prove 300, Spectroquant Prove 600		
	Parameters:	Rated voltage: 24 VDC Rated frequency: 2.5 A		
	Tested according to:	UL 61010-1:2012/R:2019-07 CSA C22.2 No. 61010-1:2012/A1:2018-11		
	This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL for USA and a Standards Council of Canada ISO/IEC 17065 accredited Certification body for Canada. Test report no.: 713298142			
Date:	2023-05-08			
				
TÜV® TÜV SÜD America, Inc. • 401 Edgewater Place Suite #600 • Wakefield • MA 01880 • USA				

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